

## Overview of all Walser® Dental Matrices:

No	Shape	Special features	Strip height	Application
1	X	also for milk molars	5,0 mm	PREMOLAR
2	X		6,5 mm	PREMOLAR
3	X		5,0 mm	MOLAR
4	X		6,5 mm	MOLAR
5	O	also for milk molars	5,0 mm	PREMOLAR-SMALL
6	O		6,5 mm	PREMOLAR-SMALL
<b>6c NEW</b>	<b>O</b>	<b>for one-sided deep caries one strip side 3 mm longer</b>	<b>5/8 mm</b>	<b>PREMOLAR-SMALL</b>
7	O	also for milk molars	5,0 mm	PREMOLAR
8	O		6,5 mm	PREMOLAR
<b>8a NEW</b>	<b>O</b>	<b>for 4-surface fillings</b>	<b>6,5 mm</b>	<b>PREMOLAR</b>
9	O		5,0 mm	MOLAR
10	O		6,5 mm	MOLAR
<b>10a NEW</b>	<b>O</b>	<b>for 4-surface fillings</b>	<b>6,5 mm</b>	<b>MOLAR</b>
<b>10c NEW</b>	<b>O</b>	<b>for one-sided deep caries one strip side 3 mm longer</b>	<b>5/8 mm</b>	<b>MOLAR</b>
<b>10d NEW</b>	<b>O</b>	<b>for one-sided deep caries one strip side 5 mm longer</b>	<b>5/10 mm</b>	<b>MOLAR</b>
11	XF	strips slit towards the end	6,5 mm	11, 12, 21, 22
12	XF	strips slit towards the end	8,0 mm	11, 12, 21, 22
13	X	extra short strips/also for milk molars	5,0 mm	13, 14, 23, 24, 31, 32, 41, 42
14	X	extra short strips	6,5 mm	13, 14, 23, 24, 31, 32, 41, 42
15	O	"	5,0 mm	31, 32, 41, 42
16	O	"	6,5 mm	31, 32, 41, 42
17	X	extra tall strips	8,0 mm	4 - 5
18	X	"	8,0 mm	6 - 8
19	O	"	8,0 mm	4 - 5
20	O	"	8,0 mm	6 - 8
21	OF	strips slit towards the end	8,0 mm	11, 21
22	OF	"	8,0 mm	12, 22
23	ON	longer strips, wide span	5,0 mm	6 - 8
<b>24</b>	<b>ON</b>	<b>attaches over cofferdam clamp!!</b>	<b>6,5 mm</b>	<b>5 - 8</b>
25	O	(Special form for milk molars; also very suitable nos. 1, 5, 7, 13)	5,0 mm	MILK MOLAR

### Sets of Walser® matrices:

**25 matrices** on sterilisation tray with/without forceps

**18 matrices** on sterilisation tray with/without forceps

**10 matrices** on sterilisation tray with/without forceps

**Kids Kit** with forceps

### Refill packs of Walser® matrices:

**Refill** of 5 matrices No. 1-22 + 25

**Refill** of 2 matrices No. 23/24

**Walser® matrix forceps** and **sterilisation trays** are available separately, too.



**Tens of thousands of dentists in more than 120 countries worldwide work with WALSER® Matrices and that number increases every day...**